

Certificate of Analysis

Redox Standard

200 ± 5mV @25°C (Pt v Ag-AgCl Ref 3M KCl)

Product No: RS200
Lot No: RS20020L1
Expiry date: 28/11/2022

Redox Potential: 205 mV @ 25°C
Date of measurement: 03/12/2020

Measurement Limits:
195 - 205mV @ 25°C

Measurement Cell:
Pt/Sample/3MKCl/AgCl/Ag

Method:
Direct potentiometry
Temperature versus redox value
T°C mV
10 225
15 216
20 208
25 200
30 191

Reference: Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 03/12/2020

Orla Sheehan

QC Technician Orla Sheehan

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-20L24