

CERTIFICATE OF ANALYSIS

Product : **HYDROGEN PEROXIDE** Analpure® **ULTRA**
 Code : UCH02161a
 Batch number : 17219071
 Expiry date : 04/02/2023 (dd/mm/yyyy)

Assay : 31 %
 Trace impurities (ppt) :

Analyte	Actual value	Analyte	Actual value
Aluminium	< 20	Neodymium	< 1
Antimony	< 1	Nickel	< 5
Arsenic	< 10	Niobium	< 0.5
Barium	< 1	Palladium	< 0.5
Beryllium	< 5	Platinum	< 1
Bismuth	< 1	Potassium	< 5
Boron	< 10	Praseodymium	< 0.1
Cadmium	< 1	Rhenium	< 0.5
Calcium	< 50	Rhodium	< 0.5
Cerium	< 0.5	Rubidium	< 1
Cesium	< 0.5	Ruthenium	< 1
Chromium	< 5	Samarium	< 0.5
Cobalt	< 1	Scandium	< 1
Copper	< 1	Selenium	< 5
Dysprosium	< 0.5	Silver	< 0.5
Erbium	< 0.5	Sodium	< 5
Europium	< 0.5	Strontium	< 1
Gadolinium	< 0.5	Tantalum	< 0.1
Gallium	< 5	Tellurium	< 1
Germanium	< 10	Terbium	< 0.1
Gold	< 5	Thallium	< 0.1
Hafnium	< 0.1	Thorium	< 0.1
Holmium	< 0.1	Thulium	< 0.1
Indium	< 1	Tin	< 5
Iron	< 20	Titanium	< 5
Lanthanum	< 0.2	Tungsten	< 1
Lead	< 0.5	Uranium	< 0.1
Lithium	< 1	Vanadium	< 1
Lutetium	< 0.1	Ytterbium	< 0.1
Magnesium	< 10	Yttrium	< 0.1
Manganese	< 1	Zinc	< 10
Mercury	< 20	Zirconium	< 0.1
Molybdenum	< 1		

Store at room temperature (15-25°C). Minimize exposure to light.



Release date: 04th February 2020

Person-in-Charge: **Vladimir Patera**

CERTIFIKÁT O ANALÝZE

Produkt : **PEROXID VODÍKU** Analpure® **ULTRA**
Objednací číslo : UCH02161a
Číslo šarže : 17219071
Exspirace : 04/02/2023 (dd/mm/rrrr)

Koncentrace : 31 %
Stopové nečistoty (ppt) :

Analyt	Aktuální hodnota	Analyt	Aktuální hodnota
Aluminium	< 20	Neodymium	< 1
Antimony	< 1	Nickel	< 5
Arsenic	< 10	Niobium	< 0.5
Barium	< 1	Palladium	< 0.5
Beryllium	< 5	Platinum	< 1
Bismuth	< 1	Potassium	< 5
Boron	< 10	Praseodymium	< 0.1
Cadmium	< 1	Rhenium	< 0.5
Calcium	< 50	Rhodium	< 0.5
Cerium	< 0.5	Rubidium	< 1
Cesium	< 0.5	Ruthenium	< 1
Chromium	< 5	Samarium	< 0.5
Cobalt	< 1	Scandium	< 1
Copper	< 1	Selenium	< 5
Dysprosium	< 0.5	Silver	< 0.5
Erbium	< 0.5	Sodium	< 5
Europium	< 0.5	Strontium	< 1
Gadolinium	< 0.5	Tantalum	< 0.1
Gallium	< 5	Tellurium	< 1
Germanium	< 10	Terbium	< 0.1
Gold	< 5	Thallium	< 0.1
Hafnium	< 0.1	Thorium	< 0.1
Holmium	< 0.1	Thulium	< 0.1
Indium	< 1	Tin	< 5
Iron	< 20	Titanium	< 5
Lanthanum	< 0.2	Tungsten	< 1
Lead	< 0.5	Uranium	< 0.1
Lithium	< 1	Vanadium	< 1
Lutetium	< 0.1	Ytterbium	< 0.1
Magnesium	< 10	Yttrium	< 0.1
Manganese	< 1	Zinc	< 10
Mercury	< 20	Zirconium	< 0.1
Molybdenum	< 1		

Skladujte při pokojové teplotě (15-25°C). Minimalizujte vystavení světlu.



Datum vydání: 04. února 2020

Odpovědná osoba: Vladimír Patera