

Certificate of Analysis

Redox Standard

358± 5mV @25°C (Pt v Ag-AgCl Ref 3M KCl)

Product No:	RS358
Lot No:	RS35821F1
Expiry date:	28/06/2023

Redox potential:	359 mV @ 25°C
Date of measurement:	03/06/2021

Measurement Limits:
353 - 363mV @ 25°C

Measurement Cell:
Pt/Sample/3MKCl/AgCl/Ag

Method:
Direct potentiometry
Temperature versus redox value

T°C	mV
10	377
15	371
20	364
25	358
30	351

Reference: Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 03/06/2021



QC Technician Orla Sheehan

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21E19