

Certificate of Analysis

Chemical Oxygen Demand - COD Reagents

COD Calibration Standard 50ppm

<i>Product No:</i>	COD50
<i>Lot No:</i>	COD5021E1
<i>Expiry Date:</i>	28/05/2023

<i>Mean Assay:</i>	50 ppm as COD
<i>Date of Measurement:</i>	31/05/2021

Specification:
49.75-50.25ppm

Method:
Standardised in accordance to in-house method, acidimetric.

Traceability:
This volumetric solution was checked by means of Sodium Hydroxide Analytical Volumetric Standard. This volumetric standard is directly traceable to a Standard Reference Material of National Institute of Standards and Technology, NIST SRM 84L, Potassium Hydrogen Phthalate.

Reference:
Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of Issue of the certificate: 31/05/2021



QC Technician Orla Sheehan

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21E16