

Certificate of Analysis

Chemical Oxygen Demand - COD Reagents

COD Calibration Standard 2000mg/l

Product No: COD2000
Lot No: COD2M20K1
Expiry date: 28/10/2022

Mean assay: 1997 mg/l as COD
Date of measurement: 14/10/2020

Specification:
1990 - 2010mg/l

Method:

Standardised in accordance with in-house method acidimetric.

Traceability:

This standard was checked by means of Sodium Hydroxide Analytical Volumetric Standard. This volumetric standard is directly traceable to a Standard Reference Material of National Institute of Standards and Technology, NIST SRM 84L, Potassium Hydrogen Phthalate

Reference:

Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 14/10/2020

Rosemary Cooney.

Q.A. Officer Rosemary Cooney

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-20K10

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com