

Certificate of Analysis

Redox Standard

400± 5mV @25°C (Pt v Ag-AgCl Ref 3M KCl)

Product No:

RS400

Lot No:

RS40021D1

Expiry date:

28/04/2023

Redox potential:

398 mV @ 25°C

Date of measurement:

06/04/2021

Specification:

395-405mV @25°C

Measurement Cell:

Pt/Sample/3MKCl/AgCl/Ag

Method:

Direct potentiometry

Temperature versus redox value

T°C mV

10 383

15 389

20 395

25 400

30 406

Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate:

06/04/2021

Eileen Malone

QC Technician Eileen Malone

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21C64

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com